



FIELD CHANGE ORDER

Service

Issued by : SSR Germany

Released : L.Küppens Groot gezeichnet

Publication No. : 4512 980 57841

Reference No. : 00 135 010

Date : April 1999

Product Group : 742



OPTIMUS RAD/RF

APPLIES TO:

All OPTIMUS RAD/RF with Rotor Control HS
9890 000 0221X 4512 104 7140X
 4512 104 3379X

TITLE:

Einschaltstrombegrenzung ROCO

LIST OF PAGES & DRAWINGS:

1 ... 2 99.0
ANR 99.0

CONTENTS:

1x power on current limitation kit 4512 104 91281

INTRODUCTION:

Symptom : Ausfälle der ROCO Einheiten
Cause :
Remedy : Einbau der Einschaltstrombegrenzung

Vorbereitungen:

OPTIMUS RAD

Vor der Änderung ist die Anlage durch Auslösen des EN-F1 spannungsfrei zu schalten und auf Spannungsfreiheit zu prüfen.

OPTIMUS RF

Vor den Serviceaktivitäten ist die Anlage spannungsfrei zu schalten. Das Auslösen des EN-F1 reicht in diesem Fall nicht aus.

Kurzbeschreibung:

Ziel dieser FCO ist es den Einschaltstrom zu begrenzen und dadurch das Anlaufgerät vor einer Zerstörung zu schützen.

MANPOWER / TIME TO COMPLETE:

1Stunde 1 Techniker

TOOLS & TEST EQUIPMENT:

Standard Werkzeugsatz

MODIFICATION KIT / PARTS REQUIRED:

Für alle o.a. Generatoren kann das folgende Kit bestellt werden:

1x power-on current limitation 4512 104 91281

Note

The modification kit will be available Free Of Charge till the end of December 2000 at the spare part centre in Hamburg.

DOCUMENTATION:

Bitte die FCO und das Manual im OPTIMUS Ordner Kapitel 8 abheften.

For local SSD use only; do not return to PMG.

4512 104 3379X

--

Not completed because customer has unit in storage.
Required parts & instructions received by the customer.

DATE _____

MAIL TO : SSD Customer Support Manager